

Integra Cloud Fax® | A2E

CONFIGURATION OPTIONS

The **Integra Cloud Fax (ICF) A2E** device fits seamlessly into a healthcare facility's existing fax infrastructure. Before install, an A2E will be enabled with the desired fax number. On location, simply plug the phone line and internet network line into the A2E device, then continue using the fax machine as normal. The A2E will then receive and send all faxes.

Faxes outgoing from the A2E will be delivered to the pharmacy through the cloud. Because that analog fax number is now connected to your Secure Exchange Network (SEN), the pharmacy is guaranteed to receive all faxes sent from that number. All faxes going from the pharmacy to the A2E number will queue in the cloud until the fax machine is ready to receive them through the A2E.

CONFIGURATION OPTIONS

Current Arrangement


Phone line is connected to a fax machine that sends and receives all faxes via the phone line. Redundant phone lines for fail over or high volume.



A2E - Same Room Configuration

Fax machine and Ethernet cords plug into the A2E device which simultaneously sends and receives faxes via the internet.*

*One A2E device is needed per fax machine.



A2E - Communications Room Configuration

Fax machine plugs into phone jack which is connected to a communications room. At that location, the phone line plugs into the A2E device along with the Ethernet cord. A2E sends and receives faxes through the cloud via the internet.

